

Saltra

MASSACHUSETTS INSTITUTE OF TECHNOLOGY  
LINCOLN LABORATORY

TO: Those Interested

FROM: J. Winett

SUBJECT: 2741 Configuration

July 10, 1969

The possible configuration for 2741 consoles is as shown below. RPQ numbers are unique to each installation even though the item is the same.

<u>Type/Feature Code</u>	<u>Name</u>	<u>Mo.</u>	<u>Rent</u>
2741-1	Communication Terminal	\$	98
*3255	Dial Up		3
*4708	Interrupt		2.50
8341	Typamatic Keys		5
9104	10 Characters/Inch Spacing		NC
*9114	103A Data Set Attach		NC
9162	13 5/8 Forms width 6 lines/inch		NC
9435	3 or 6 lines/Inch		NC
9509	Pinfeed Platen (SUC \$59.25)		
9571	Type Element PTTC/EBCD		NC
9881	115 Volt AC 60 Cycle		NC
+ ** *E46766	Receive Int. Control (RPQ)		8
+ ** E46151	Print Inhibit (RPQ)		10
+ ** E67038	Ribbon Shift		10
+ **E67040	Auto address answerback		12
+ **	Auto EDT	Total	\$148.50

\*Required for CP/CMS

\*\*Required for CTSS

+ Required for Multics

J. Winett

JMW:lg

FOR LABORATORY USE ONLY